



Cellevate

Immunocytochemistry of cells in nanofiber scaffolds.*

1. Wash cells or spheroids cultured in the nanofiber scaffolds using preheated PBS (with 1 mM CaCl₂ and 1 mM MgCl₂).
2. Fix cells using a 4 % paraformaldehyde in PBS (pH 7.4) for 10 min at room temperature. 0.1 % glutaraldehyde may be added to enhance cell fixation.
Note: Avoid using organic solvents that may deteriorate the nanofiber scaffold, e.g. methanol or acetone.
3. Rinse cultures 3 x 5 min with PBS.
4. Cells may be permeabilized and non-specific binding blocked by incubating cultures for 1 h at room temperature using 0.3 % (v/v) Triton X-100 and 5 % (w/v) BSA (or other suitable blocking agent) in 100 mM PBS.
Note: Increase blocking agent concentration if trouble with the nanofibers adsorbing the antibodies occur.
5. Rinse cultures 3 x with PBS.
6. ICC can now be performed using standard protocols.
Note: For high-magnification microscopy (>10X) see the protocol "*Mounting nanofiber scaffolds for microscopy*" for suggestions on how to mount your samples.

* Suggested procedure, please adjust according to your experimental needs.