

## Immunocytochemistry of cells in nanofiber scaffolds\*

1. Wash cells or spheroids cultured in the nanofiber scaffolds using preheated PBS (with 1 mM CaCl<sub>2</sub> and 1 mM MgCl<sub>2</sub>).
2. Fix cells using a 4 % paraformaldehyde in PBS (pH 7.4) for 10 min at room temperature. 0.1 % glutaraldehyde may be added to enhance cell fixation.  
**Note:** Avoid using organic solvents that may deteriorate the nanofiber scaffold, e.g. methanol or acetone.
3. Rinse cultures 3 x 5 min with PBS.
4. Cells may be permeabilized and non-specific binding blocked by incubating cultures for 1 h at room temperature using 0.3 % (v/v) Triton X-100 and 5 % (w/v) BSA (or other suitable blocking agent) in 100 mM PBS.  
**Note:** Increase blocking agent concentration if trouble with the nanofibers adsorbing the antibodies occur.
5. Rinse cultures 3 x with PBS.
6. ICC can now be performed using standard protocols.  
**Note:** For high-magnification microscopy (>10X) see the protocol "Mounting nanofiber scaffolds for microscopy" for suggestions on how to mount your samples.

\*Suggested procedure, please adjust according to your experimental needs.