



Cellevate

Cell recovery by trypsinisation.*

Cells can be extracted from the scaffold by trypsinisation, using accutase, or other appropriate enzyme, and gentle agitation using standard protocols. Avoid using cell scrapers. Instead, gently dislodge the cells by careful pipetting and tapping on the culture plates. Additionally, in general, longer incubation times and more rigorous shaking is required compared to a 2D culture.

Please optimise your trypsinisation protocol for each cell type to achieve full cell recovery.

* Suggested procedure, please adjust according to your experimental needs.